

ZUNTO

Thin and light, the IM02 medical pen lamp is designed for night surveillance in bedrooms and assessment of pupillary reflexes

Discreet

Its discreet beam that respects patients' sleep is ideal for visiting the night room

Pen Format

Its pen format allows hanging on any type of pocket or in your first aid kit



Luminous flow for pupillary osculation



Operating time



Lighting distance



Designer and manufacturer, Lagolight offers solutions of portable lighting for professionals

Not just a lamp, this is a tool

ZUNTO



Leclanché becomes Zunto

IM02 medical pen lamp

Reference: YL000IM023

Lighting

- 1 incandescent bulb with optical bulb
- Pen type push button ignition
 - Lighting intensity: 10 lumens.
 - Lighting distance: 10 meters.

Energy – Operating time

- 2 LR3/AAA alkaline 1.5V batteries (not supplied).
- Operating time: 13 hours.
- Compatible with rechargeable RC3/AAA 1.2V accumulators (not supplied).

Application – Use

- Rescue - Oscultation - Medical professions - auxiliary lamp ...
- Operating temperature: -25°C / +50°C.

Lamp enclosure

- Lacquered aluminium case.
- Stainless steel contact spring.
- Lamp cleaning by rinsing with clear water and drying.

Weight – Dimensions

- Extremely light weight: 32 grams (batteries included).
- Dimensions: length 140 mm x diameter 12.4 mm.
- Pen hook for hanging on a pocket.
- Case color: white with a black cap.

Warranty

- Warranty: 2 years.



Assessment of pupillary reflexes



LAGOLIGHT
Headlamps & Portable lighting

531 route des Vernes 74370 Pringy - France
+33 (0)4 50 60 55 92 contact@lagolight.com
www.lagolight.com

Designer &
manufacturer
for
30 years